



Probe Introducer – 10G and 12G
(ECCOAXCAN-10 & 12G)








Vacuum Package
(SCVACPAK3000)



Rinse Tubing Set for Auto Rinse Option
(RINSETUB)




Sample Sizes



Probe Size	7G	10G	12G
Probe Length	34.3 cm	34.2 cm	34.1 cm
Needle Length	11.4 cm	11.3 cm	11.2 cm
Needle Diameter	5.0 mm	4.0 mm	3.4 mm
Dead Space	10.8 mm	9.9 mm	8.7 mm
Sample Notch	19.0 mm	19.0 mm	19.0 mm
Tip to Center Notch	20.3 mm	19.4 mm	18.2 mm
Fill Volume	3.5 cc	2.0 cc	1.2 cc


EnCor[®]
360



Tri-Concave Trocar Tip



Coaxial Cannula with 1cm depth markings



360° Sample Chamber



Sample Collection



Vacuum Package
(SCVACPAK3000)



3 Morgan
Irvine, California 92618
877-210-1300 Customer Service
949-362-3519 Fax
customerservice@senorx.com
www.senorx.com

Rx Only
Please consult product labels and inserts for any indications, contraindications, hazards, warnings, cautions and directions for use.

©Copyright 2009 Senorx Inc. All rights reserved.
EnCor is a registered trademark of Senorx Inc.