

## Proven Performance For You. Enhanced Comfort For Her.

- ✓ Vacuum ports enable the removal of fluid and air and facilitate tissue conformance to the balloon for improved tissue dosing
- ✓ Multiple offset lumens optimize dosimetric coverage of targeted tissue while minimizing radiation dose to skin and ribs
- ✓ Polyurethane balloon maintains symmetry and reduces potential for balloon rupture

Inflation Luer

Vacuum Luer

Stiffening Stylet

Radiation Source Lumens

5 Treatment Lumens

Vacuum Ports

Orientation Line

Introducing ConturaFlex. Now our proven MLB design is more flexible for even greater APBI patient comfort.

**NEW!** 50% more flexible shaft ensures greater patient comfort

**NEW!** Shorter distal tip reduces spacing for enhanced tissue conformance and treatment planning

# Contura<sup>®</sup> Flex

50% More Flexible Shaft  
For Enhanced Patient Comfort



## ConturaFlex MLB Ordering Information

Product Number	Product Description
<b>ConturaFlex MLB – Multi-Lumen Balloon</b>	
<b>B112-45</b>	<b>4.0-5.0 cm</b> ConturaFlex MLB (B013-45) with <b>Tri-Concave™ Tip</b> Introducer (B900-10) and Procedure Pak
<b>B112-56</b>	<b>4.5-6.0 cm</b> ConturaFlex MLB (B013-56) with <b>Tri-Concave™ Tip</b> Introducer (B900-10) and Procedure Pak
<b>B122-45</b>	<b>4.0-5.0 cm</b> ConturaFlex MLB (B013-45) with <b>Pyramidal Tip</b> Introducer, (B900-17) and Procedure Pak
<b>B122-56</b>	<b>4.5-6.0 cm</b> ConturaFlex MLB (B013-56) with <b>Pyramidal Tip</b> Introducer, (B900-17) and Procedure Pak
<b>B132-45</b>	<b>4.0-5.0 cm</b> ConturaFlex MLB (B013-45) with <b>Blunt Tip</b> Introducer (B900-09) and Procedure Pak
<b>B132-56</b>	<b>4.5-6.0 cm</b> ConturaFlex MLB (B013-56) with <b>Blunt Tip</b> Introducer (B900-09) and Procedure Pak
<b>B113-45</b>	<b>4.0-5.0 cm</b> ConturaFlex MLB (B013-45) and Procedure Pak ( <b>No Introducer</b> or split sheath)
<b>B113-56</b>	<b>4.5-6.0 cm</b> ConturaFlex MLB (B013-56) and Procedure Pak ( <b>No Introducer</b> or split sheath)
<b>For reference only, not ordered separately.</b>	
<b>4.0 – 5.0 cm Procedure Pak</b> includes drainage catheter, <b>three</b> 30 ml inflation syringes, <b>one</b> 10 ml inflation syringe, <b>one</b> #11 scalpel, contrast media tray, radiation lumen caps	
<b>4.5 – 6.0 cm Procedure Pak</b> includes drainage catheter, <b>five</b> 30 ml inflation syringes, <b>one</b> 10 ml inflation syringe, <b>one</b> #11 scalpel, contrast media tray, radiation lumen caps	



INNOVATIVE BREAST CARE AT EVERY STEP

Senorx, Inc.  
3 Morgan  
Irvine, CA 92618-1917  
877-210-1300  
949-362-3519 Fax  
customerservice@senorx.com  
www.senorx.com

*Rx Only*  
Please consult product Directions For Use for indications, contraindications, warnings, and precautions.

Contura is a registered trademark of Senorx, Inc. Tri-Concave is a trademark of Senorx, Inc.

©Copyright 2010 Senorx, Inc. All rights reserved.

